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
## CONTENTS

1. Preface-----	1
2. Main technical indexes-----	2
3. Schematic drawing of the machine-----	2
4. Preparation before using the machine-----	4
5. Operating instructions-----	5
6. Maintenance-----	5
7. Wiring diagram-----	7

## SAFETY MEASURES

- » Read this manual carefully before installing and using the equipment.
- » In the meantime, follow the warnings, especially the recommended content.
- » Before using the equipment, check the input voltage range that complies with the requirements. The operator must be trained, and once qualified, you can operate this machine.
- » Check that the surrounding environment complies with the installation requirements.
- » Operate the equipment in accordance with the regulations.
- » Do not allow the equipment to be disassembled when a fault occurs.
- » You should invite a professional to repair the equipment to avoid shock.
- » When cleaning the equipment, the operator should disconnect the power supply. Keep this instruction manual.

## FOREWORD

- » The MLP-CM60 sample cutting machine is suitable for cutting a variety of samples, metallic and non-metallic materials to observe petrographic structure.
- » The machine has a cooling device. When using configured coolant, heat can be dissipated to prevent overheating.
- » The operation of the machine is easy for the user.

Model	MLP-CM60
Voltage	380V, 50Hz
Spindle Speed (r/min)	2800
Cutting wheels size (mm)	Φ250×1.5×Φ25.4
Max Cutting Capacity (mm)	Φ60
Power (KW)	2.2
Dimension (mm)	830×620×450
Weight (kg)	82

Fig.1


Fig.2


Fig.2

## PREPARATION BEFORE WORK

4.1 The machine is installed on a stable work table.

4.2 The operator first installs the supply and drainage pipes to ensure that there are no leaks. The drainage pipe should be lower than the inlet pipe to prevent poor drainage.

4.3 Check the power lines, plugs and other electrical components and then connect the power supply to start the machine. On the website The spindle direction should be in the same direction as indicated by the machine.

4.4 Loading and unloading the cutting blade

1) First, open the cover cover, press the "cutting button control" button to the position as shown in Figure 4, and then press the locking button with the right hand. Turn the spindle with your left hand until the lock button locks into place. Finally, open the left-hand nut bolt with the supplied wrench and the clamping and cutting disc.

2) Install the clamps and cutting discs into the spindle. Tighten the left-hand nut with the wrench. Then additionally the cutting discs. If the angular momentum is greater than 2 mm, the operator should check the cutting wheel.

3) After replacing the 250 mm diameter cutting wheel, be sure to adjust the "cutting stroke control handle" to the position as shown in Figure 5. The "cutting stroke control handle" can only be adjusted to the position shown in Figure 6 to increase the consumption of the blade when the consumption of the blade reaches 215 mm or when a cutting blade smaller than 215 mm is replaced.

Fig.4

Fig.5

Fig.6

- ⚠ 1. Remember to turn off the „off switch” or press the „emergency button” when changing the cutting wheel.
- 2. Damage to the machine caused by improper adjustment of the "cutting stroke control" is not covered under warranty.

4.5 Adjusting the handle and cutting stroke

» First, adjust the relative position of the cutting handle with the Cutting Handle Arresting Control to confirm the cutting comfort. The handle is controlled via the cutting stroke (see Figure 2).

4.6 Installing the workpiece

» Correct holding of the sample is an important guarantee for safe operation. The operator first inserts the workpiece into the jaws of the vise.

» He then presses the handle to touch the workpiece and secures the workpiece via the handle.

» Finally, check that the workpiece is stable.

## OPERATING INSTRUCTIONS

1. The cutting force should be uniform. Do not apply too much pressure - it is recommended to cut 1-2 mm with each stroke and lift the handle slightly.
2. Be sure to open the ball valve on the water pump before starting work. The valve should be open when cutting the workpiece; it is necessary to close the ball valve on the pump before cleaning.
3. The coolant must have good lubricity. It is recommended to use the coolant supplied with the machine (the first batch is free) or it can be purchased from our company when exhausted. Ordinary tap water should not be used as coolant.
4. Stoneware, ceramics and similar materials must be cut with diamond cutting wheels.

### 5.1 Operating keys: (see figure 1)

- 1) Cleaning: a separate "Self Cleaning" key and a "Cleaning" key are used to operate the pump. When this key is pressed, the indicator light will come on and the water pump will enter the operating state.
- 2) Starting: reset button, start button. During normal operation, this button is lit.
- 3) Stop: Press the button to stop the machine during the incomplete process. If the button is not active, the light remains off.
- 4) Emergency stop: Emergency stop button for emergency, press this button for protection. (Meanwhile, the indicator of this emergency stop button is not lit.)


### 5.2 Operating Instructions:

- 1) Complete the preparation before starting operation. Cover the machine and turn on the power switch. As soon as the control panel indicator lights up, the machine is ready for use.
- 2) Open the ball valve on the water pump and press the start button. The motor will start. The blade will touch the sample during cutting. When cutting is complete, rotate the swing arm handle back to its original position. Then press the stop button - the motor will switch off
- 3) Close the valve on the water pump, press the cleaning button and clean the cutting area.

1. The cutting force should be uniform and not too high. It is recommended to cut the material in layers 1-2 mm thick and to lift the handle slightly.
2. Before starting to cut, it is done that you open the ball valve on the waterjet. When you have finished cutting, close it before starting the cleaning process.
3. Use coolants with good lubricity (one cartridge is included with the machine and can be used by our company). Ordinary tap water is not suitable as coolant.
4. Use only diamond cutting wheels for cutting rocks, ceramics and similar materials.

## MAINTENANCE

- 6.1 The machine must be fitted with a reliable earthing line.
- 6.2 Inlet and outlet piping shall be in good condition but shall not leak.
- 6.3 Do not allow broken cutting wheels to cut the sample.
- 6.4 Do not allow the use of cutting wheels whose speed is less than 35 m/s.
- 6.5 After operating for a time, the operator must clean the machine to make the machine have a longer service life.
- 6.6 Change coolants regularly to prevent them from deteriorating.
- 6.7 When you detect abnormal sound, you should immediately stop the machine and check.
- 6.8 Replenish the lubricating oil regularly.

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**INSIZE is a global manufacturer of measuring technology** bringing technological innovations with representation in 75 countries around the world. Measuring instruments of the INSIZE brand represent the optimal solution without compromises to meet even the most demanding measuring needs.

**INSIZE measuring instruments will pleasantly surprise you with:**  
| high-quality craftsmanship | reliable performance | affordable prices

INSIZE provides a complete portfolio of measuring instruments » calipers, height gauges, dial indicators, micrometers, roughness testers, hardness testers, measuring microscopes, optical measuring devices, profile projectors, testing machines, contour gauges, roundness measuring machines, thickness gauges, gages, protractors, force gauges, meters, scales, videoscopes, torque wrenches and accessories including stands for measuring instruments.

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